

Quick Start Guide Aco-Dyes™ for Labelling Extracellular Vesicles

(A) Dye Reconstitution

① Add 100µL of filtered
 1X PBS / water



Aco-Dye
 100 Tests Size

Photophysical Properties of Aco-Dyes

Product	Aco-430™	Aco-490™	Aco-520™	Aco-600™	Aco-650™	Aco-800™
$\lambda_{abs} (SUUV)$ (nm)	369	422	465	525	632	735
$\lambda_{em} (SUUV)$ (nm)	403-460	458-508	511-530	586-635	645-661	775-818
Recommended Laser(s)	UV (355nm)	Violet (405nm)	Blue (488nm)	Blue/Yellow (488/561nm)	Red (638nm)	Red (638nm)

Ensure that your instrument is equipped with the appropriate lasers and filters to detect the chosen dye(s) before proceeding with your experiment.

Storage Conditions


- Dry Dyes: Room Temperature
Refer to expiry date on tube
- Reconstituted Dyes: 2-8°C
Up to 1 month
- Frozen Reconstituted Dyes: -80°C
Up to 3 months; avoid freeze-thaw

Aco-Dyes perform best when freshly reconstituted. Perform Steps 2 to 5 prior to using previously reconstituted dyes.


② Vortex to mix



③ Quick spin



④ **[Important!]**
 Sonicate at 40°C
 for 15 minutes



If a sonicator is not available, incubate at room temperature for 15 minutes.

⑤ Repeat Steps 2-3:
 Vortex > Quick spin




25µM dye
 stock solution


(B) EV Labelling for Flow Cytometry

**It is essential to optimise the staining concentration for your EV type, dye and instrument.*


⑥ Dilute dye stock to 10µM
 e.g. 9µL buffer + 6µL dye



⑦ Incubate 1e10 particles/mL EV sample with 1µM* final dye concentration for 1 hour at 37°C, protected from light
 e.g. 45µL of 1e10 particles/mL EV + 5µL of 10µM dye



⑧ Dilute stained EV samples just before flow cytometry acquisition to achieve instrument manufacturer's recommended event rate



Important Considerations

1. It is recommended to include the following **controls** to identify background signals and rule out false positives.
 - **Dye Only** [*same dye conc as stained EV*]
 e.g. 45µL buffer + 5µL of 10µM dye
 - **EV Only**
 e.g. 45µL of 1e10 particles/mL EV + 5µL buffer
 - **Buffer Only**
2. The **dilution factor** optimised in Step 8 for the stained EV samples should be applied uniformly to all related samples, including controls.
3. Using the **same instrument settings** (gain, threshold and width) and **acquisition time** across all samples ensures reliable comparisons.